



## Features

- Resin/ Rubber for Medical use (polycarbonate / Silicone rubber)
- Assembled in a Clean room in Japan.
- Tube connection (Recommended Inside diameter:  $\Phi 2.5\text{mm}$ )
- Compact (Outside diameter  $\Phi 11\text{mm}$  × Length: 23.1mm)

## Specifications (tentative)

- Cracking Pressure: 2kPa (after priming wash)
- Min. opening differential pressure: 15kPa (water)
- Pressure resistance: 150kPa ~
- Flow rate: 550mℓ~/min (100 kPa)

Note) Standard: No sterilization

**PATENTED**



**We can develop  
“Resin + Rubber” Components based on your request.**

## 【Evaluation method and results】

### ● Opening valve

Water and air are applied from the forward direction at 2kPa, and it is confirmed that the valve opened. (after priming wash)

### ● Check valve (JIS T 3211 5.16 Standard)

After priming, pressure from 15 kPa to 150 kPa are applied from the check direction in increments of 5 kPa, and after observed for 15 seconds, it is confirmed that there is no continuous backflow.


### ● Pressure resistance (JIS T 3248 5.5.1 Standard)

The check valves are submerged in a container filled with water at  $37\pm 1^{\circ}\text{C}$ , and air pressure of 150kPa is applied from the check direction for 10 minutes. It is confirmed that no continuous bubbles are generated.

### ● Flow rage

Random sample	1	2	3	4	5	6	7	8	9	10	Ave.
Water[ml/min]	600	610	590	600	610	595	610	580	615	590	600

\*Inside diameter  $\phi 2.5\text{mm}$  and inlet pressure 100kPa



< Contact >



**FUJIKURA COMPOSITES**

GOTANDA JP BUILDING 4F  
8-4-13 Nishigotanda, Shinagawa-ku, Tokyo, Japan 141-0031  
TEL. +81-3-5747-9223  
FAX. +81-3-5747-9217

Web: <https://www.fujikura-control.com/english>  
Founded: 1901

Sales: \$281 million US dollar (2023 March, consolidated)  
ISO certified: 9001, 14001, 13485

E-mail: [seigyotoiawase.en@fc.fujikura.co.jp](mailto:seigyotoiawase.en@fc.fujikura.co.jp)